

# T.gondii ROP4 Recombinant Protein



RPPB5544

---

## Product Information   Protein Information

**Product SKU:**

RPPB5544

**Host:**

Escherichia Coli.

**Protein description:**

The E.coli derived recombinant artificial mosaic protein contains the ROP4 (RH2) immunodominant regions, the protein is fused to Glutathione transferase.

**Formulation:**

25mM Tris-HCl pH 8.0, 5mM Glutathione, 1.5M urea and 50% glycerol.

**Purity:**

Protein is >95% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

Toxoplasma protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.