

# T.pallidum p15 (Partial) Recombinant Protein



RPPB5545

---

## Product Information    Protein Information

**Product SKU:**

RPPB5545

**Host:**

Escherichia Coli.

**Protein description:**

The E.Coli derived recombinant protein is fused at N-terminus with GST tag and contains the Trp. Pallidum p15 immunodominant regions.

**Formulation:**

70mM Tris-HCl pH-8, 84mM NaCl, 14mM Glutathione, 30% Glycerol & 0.2% Sarcosil.

**Purity:**

Treponema Pallidum protein is >95% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

Treponema Pallidum protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.