

T.pallidum p15 Recombinant Protein



RPPB5546

Product Information Protein Information

Product SKU:

RPPB5546

Host:

Escherichia Coli.

Protein description:

The E.Coli derived recombinant protein contains the Trp. Pallidum p15 immunodominant regions. The protein contains beta-galactosidase (114 kDa) fused at the N-terminus.

Formulation:

20mM Tris-HCl pH-8, 10mM B-ME and 8M urea.

Purity:

Treponema Pallidum protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

Protein should be stored at 4°C.Refrigerate Upon arrival. DO NOT FREEZE.