

# T.pallidum p15 (Partial) Recombinant Protein



RPPB5547

---

## Product Information   Protein Information

**Product SKU:**

RPPB5547

**Host:**

Escherichia Coli.

**Protein description:**

The E.Coli derived recombinant 6xHis tag fusion protein is a multimer having a molecular mass of 15kDa and contains the Trp. Pallidum p15 immunodominant regions and six histidines fused at the C- terminus.

**Formulation:**

70mM Tris-HCl pH 8.0, 50mM NaCl, 50% Glycerol, 1.5 M Urea.

**Purity:**

Treponema Pallidum protein is >95% pure as determined by SDS- PAGE.

**Stability:**

Treponema Pallidum protein should be stored at 2-8°C. Do NOT freeze.