

T.pallidum p17 (Partial) Recombinant Protein



RPPB5549

Product Information Protein Information

Product SKU:

RPPB5549

Host:

Escherichia Coli.

Protein description:

The E.Coli derived recombinant protein is fused at N-terminus with 6xHis tag and contains the Trp. Pallidum p17 immunodominant regions.

Formulation:

70mM Tris-HCl pH8.0, 50mM NaCl, 50% Glycerol, 1.5M Urea.

Purity:

Treponema Pallidum protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

Treponema Pallidum protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.