

# T.pallidum p41 Recombinant Protein



RPPB5552

---

## Product Information    Protein Information

**Product SKU:**

RPPB5552

**Host:**

Escherichia Coli.

**Protein description:**

The E.coli derived recombinant protein contains the outer membrane T.Pallidum p41 immunodominant regions. The protein is fused with a Beta-galactosidase tag.

**Formulation:**

10mM Tris-HCl pH8.0, 1mM EDTA, 1mM DTT and 8M urea.

**Purity:**

Treponema pallidum protein is >90% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

T.pallidum p41 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.