

T.pallidum p41 Mosaic Recombinant Protein



RPPB5553

Product Information Protein Information

Product SKU:

RPPB5553

Host:

Escherichia Coli.

Protein description:

The E.coli derived recombinant protein contains the outer membrane T.Pallidum p41 immunodominant regions. The protein is fused with a GST tag.

Formulation:

25mM Tris-HCl pH-8, 60mM NaCl & 50% glycerol.

Purity:

Treponema Pallidum protein is >90% pure as determined by 10% PAGE (coomassie staining).

Stability:

Treponema Pallidum protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.