

T.pallidum p47 Recombinant Protein



RPPB5554

Product Information Protein Information

Product SKU:

RPPB5554

Host:

Escherichia Coli.

Protein description:

The E.Coli derived recombinant protein contains the T.Pallidum p47 immunodominant regions. The protein contains beta-galactosidase (114 kDa) fused at the N- terminus.

Formulation:

8M urea, 10mM Tris-HCl pH 8.0, 1mM EDTA and 1mM DTT.

Purity:

Treponema Pallidum protein is >95% pure as determined by SDS-PAGE.

Stability:

Treponema Pallidum although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.