

T.pallidum p47 (Partial) Recombinant Protein



RPPB5556

Product Information Protein Information

Product SKU:

RPPB5556

Host:

Escherichia Coli.

Protein description:

The E.Coli derived recombinant protein is fused at the C- terminus with a 6xHis Tag and contains the T.Pallidum p47 immunodominant regions.

Formulation:

70mM Tris-HCl pH8.0, 50mM NaCl, 50% Glycerol and 1.5M urea.

Purity:

Treponema Pallidum protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

Treponema Pallidum although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.