

T.pallidum TmpA partial Recombinant Protein



RPPB5557

Product Information Protein Information

Product SKU:

RPPB5557

Host:

Escherichia Coli.

Protein description:

The E.Coli derived 68kDa recombinant protein contains the Trp. Pallidum TmpA immunodominant regions, 23-41 amino acids and 288-325 amino acids.

Formulation:

10mM Tris-HCl pH8.0, 20mM DDT, 1mM EDTA and 8M urea.

Purity:

Protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

T.pallidum TmpA partial although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.