

T.pallidum TmpA Recombinant Protein



RPPB5558

Product Information Protein Information

Product SKU:

RPPB5558

Host:

Escherichia Coli.

Protein description:

The E.Coli derived recombinant protein contains the T.Pallidum TmpA immunodominant regions. The protein contains beta- galactosidase (114 kDa) fused at the N- terminus.

Formulation:

8M urea, 20mM Tris-HCl pH-8 & 10mM B-ME.

Purity:

Treponema Pallidum protein is >90% pure as determined by 10% PAGE (coomassie staining).

Stability:

Treponema Pallidum protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.