

Treponema TmpA 45kDa Recombinant Protein



RPPB5559

Product Information Protein Information

Product SKU:

RPPB5559

Host:

Escherichia Coli.

Protein description:

Recombinant Treponema TmpA produced in E.coli is a non-glycosylated polypeptide chain having a molecular mass of 45kDa and fused to a His tag at N-terminus.

Appearance:

Sterile Filtered White lyophilized (freeze-dried) powder.

Formulation:

Lyophilized from 1mg/ml in 20mM sodium carbonate pH-10.

Purity:

Greater than 95.0% as determined by SDS-PAGE.

Solubility:

It is recommended to reconstitute the lyophilized Treponema TmpA in sterile 18M-cm H₂O not less than 100µg/ml, which can then be further diluted to other aqueous solutions.

Stability:

Treponema TmpA although stable at room temperature for 4 weeks, should be stored below -18°C. Please prevent freeze thaw cycles.