

# HIV-1 p24 Core Recombinant Protein



RPPB5563

---

## Product Information    Protein Information

**Product SKU:**

RPPB5563

**Host:**

Escherichia Coli.

**Protein description:**

HIV-1 p24 Core is a 24kDa non-glycosylated polypeptide chain, containing the full length sequence of the HIV-1 p24.

**Appearance:**

Sterile filtered colorless clear solution.

**Formulation:**

50mM Tris, pH 8.5 and 8M urea.

**Purity:**

Greater than 95.0% as determined by HPLC analysis and SDS-PAGE.

**Stability:**

HIV-1 p24 core although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.