

# HIV-1 p24 Recombinant Protein



RPPB5565

---

## Product Information   Protein Information

**Product SKU:**

RPPB5565

**Host:**

Escherichia Coli.

**Protein description:**

The E.coli derived 39 kDa biotin labeled recombinant protein is a non-glycosylated polypeptide chain, containing the HIV-1 p24 immunodominant regions, amino acids 77-436. The HIV-1 p24 Biotinylated is fused to beta-galactosidase (114 kDa) at the N-terminus.

**Appearance:**

Sterile filtered colorless clear solution.

**Formulation:**

8M Urea, 20mM Tris-HCl pH 8.0, 10mM b-mercaptoethanol.

**Purity:**

Greater than 95.0% as determined by SDS-PAGE.

**Stability:**

HIV-1 p24 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.