

HIV-1 p24 Recombinant Protein



RPPB5566

Product Information Protein Information

Product SKU:

RPPB5566

Host:

Escherichia Coli.

Protein description:

The E.coli derived HRP recombinant protein is a non-glycosylated polypeptide chain, containing the HIV-1 p24 immunodominant regions and fused to a GST tag at N-terminus.

Appearance:

Sterile filtered colorless clear solution.

Formulation:

1.5M urea, 25mM Tris-HCl pH-8.0 and 50% Glycerol.

Purity:

Greater than 90.0% as determined by HPLC analysis and SDS-PAGE.

Stability:

HIV-1 p24 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.