

HIV-1 p24 Recombinant Protein



RPPB5568

Product Information Protein Information

Product SKU:

RPPB5568

Protein description:

HIV-1 p24 recombinant is a full-length recombinant protein having a molecular mass of 24 kDa.

Host:

Escherichia Coli.

Appearance:

Sterile filtered white lyophilized powder.

Formulation:

HIV-1 p24 was lyophilized from a solution containing 1% glycerol.

Purity:

Greater than 95.0% as determined by SDS-PAGE.

Solubility:

Reconstitute with aprotogenic sterile water or 1X PBS.

Stability:

HIV-1 p24, 24 kDa although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.