

# HIV-1 gp41/120 Recombinant Protein



RPPB5577

---

## Product Information    Protein Information

**Product SKU:**

RPPB5577

**Protein description:**

HIV-1 gp41-gp120 is a non-glycosylated polypeptide chain fused to 6xHis Tag.

**Host:**

Escherichia Coli.

**Appearance:**

Sterile filtered colorless clear solution.

**Formulation:**

150mM NaCl, 25mM sodium phosphate pH-8, 50% glycerol and 150mM imidazole.

**Purity:**

Greater than 95.0% as determined by SDS-PAGE.

**Stability:**

HIV-1 gp41/120 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.