

HIV-1 gag p17 Recombinant Protein



RPPB5582

Product Information Protein Information

Product SKU:

RPPB5582

Host:

Escherichia Coli.

Protein description:

HIV-1 p17, p24 is a non-glycosylated polypeptide chain, containing sequence of HIV-1 immunodominant regions p17-p24. The protein is fused to a GST tag at N-terminus.

Appearance:

Sterile filtered colorless clear solution.

Formulation:

1.5M urea, 25mM Tris-HCl pH 8.0, 0.2% Triton-X & 50% Glycerol.

Purity:

Greater than 95.0% as determined by SDS-PAGE.

Stability:

HIV-1 gag p17, p24 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.