

HIV-1 gag p17 Recombinant Protein



RPPB5583

Product Information Protein Information

Product SKU:

RPPB5583

Host:

Escherichia Coli.

Protein description:

HIV-1 p17, p24, gp120 is a 70 kDa non-glycosylated polypeptide chain, containing sequence of HIV-1 immunodominant regions p17,p24, gp120, the protein is fused to GST at N-terminus.

Appearance:

Sterile filtered colorless clear solution.

Formulation:

8M urea, 25mM Tris-HCl pH 8.0, 5mM b-mercaptoethanol.

Purity:

Greater than 95.0% as determined by HPLC analysis and SDS-PAGE.

Stability:

HIV-1 gag p17,p24, gp120 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.