

HIV-1 gag p17-p24 Recombinant Protein



RPPB5584

Product Information Protein Information

Product SKU:

RPPB5584

Host:

Escherichia Coli.

Protein description:

HIV-1 -p17-p24, gp-41-gp120 is a non-glycosylated polypeptide chain, containing sequence of HIV-1 immunodominant regions p17-p24, gp41-gp120.

Appearance:

Sterile filtered colorless clear solution.

Formulation:

20mM PBS pH7.8, 20mM NaCl, 1mM DTT & 8M urea.

Purity:

Greater than 95.0% as determined by HPLC analysis and SDS-PAGE.

Stability:

HIV-1 gag p17-p24, gp41-gp120 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.