

HSV 2 gG Recombinant Protein



RPPB5594

Product Information Protein Information

Product SKU:

RPPB5594

Protein description:

The HSV-2 gG protein is a synthetic protein which containing the HSV-2 gG immunodominant regions.

Appearance:

Sterile Filtered clear solution.

Formulation:

HSV-2 gG solution contains 25mM Tris-HCl pH8.0.

Purity:

Protein is >99% pure as determined by Analytical HPLC and Mass Spec.

Stability:

HSV-2 gG protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.