

HIV-2 gp32 Recombinant Protein



RPPB5598

Product Information Protein Information

Product SKU:

RPPB5598

Host:

Escherichia Coli.

Protein description:

HIV-2 gp32 Biotin Labeled recombinant- contains the full-length sequence of HIV-2 envelope immunodominant regions gp32 having a Mw of 32kDa and fused to a beta-galactosidase at N-terminus.

Appearance:

Sterile filtered colorless clear solution.

Formulation:

0.01M Na₂CO₃, 10mM EDTA, 14mM beta-ME and 0.02% Sarcosyl.

Purity:

Greater than 95.0% as determined by SDS-PAGE.

Stability:

HIV-2 gp-32 although stable at room temperature for 3 weeks, should be stored at 4°C.