

HTLV-1 p24 core Recombinant Protein



RPPB5606

Product Information Protein Information

Product SKU:

RPPB5606

Protein description:

The E.Coli derived recombinant protein contains the full-length sequence of HTLV-I p24, spanning all of p24. The protein contains 214 amino acids having an Mw of 24kDa.

Formulation:

50mM NaPO₄ pH 6.0, containing 1mM DTT & 1mM EDTA.

Purity:

HTLV-1 p24 protein is >95% pure as determined by 10% PAGE (coomassie staining) and RP-HPLC.

Stability:

HTLV-1 p24 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.