

# HTLV-1 p24 Recombinant Protein



RPPB5607

---

## Product Information   Protein Information

**Product SKU:**

RPPB5607

**Host:**

Escherichia Coli.

**Protein description:**

HTLV-1 p24 Recombinant produced in E.Coli covers full length of HTLV-1 p24 and contains 188 amino acids having a molecular size 21kDa. HTLV-1 p24 is purified by proprietary chromatographic technique.

**Appearance:**

Sterile Filtered clear solution.

**Formulation:**

HTLV-1 p24 protein solution containing PBS, pH-7.4.

**Purity:**

Protein is >95% pure as determined by 12% SDS-PAGE (coomassie staining).

**Stability:**

HTLV-1 p24 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.