

HTLV-1 gp21 Recombinant Protein



RPPB5608

Product Information Protein Information

Product SKU:

RPPB5608

Protein description:

The E.Coli derived recombinant protein contains the gp21 immunodominant region, 351-404 amino acids, the MW is 21kDa. Fused with a GST tag at N-Terminus.

Formulation:

60mM NaCl, 50mM Tris-HCl, pH-8.0, 0.25% Sarkosil, 50% Glycerol and 10mM Glutathione.

Purity:

HTLV-1 gp21 is >95% pure as determined by 10% SDS-PAGE (coomassie staining) and HPLC-C4.

Stability:

HTLV-1 gp21 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.