

HTLV-1 mosaic Recombinant Protein



RPPB5610

Product Information Protein Information

Product SKU:

RPPB5610

Protein description:

The E.Coli derived recombinant mosaic protein contains the gp21 and gp46 immunodominant regions, 374-400 amino acids and 190-207 amino acids, the MW is 42kDa. The protein is fused with GST at N-terminus.

Formulation:

60mM NaCl, 50mM Tris-HCl pH 8.0, 0.25% Sarkosil, 50% glycerol and 10mM glutathione.

Purity:

HTLV-1 Mosaic protein is >95% pure as determined by 10% PAGE (coomassie staining) and RP-HPLC.

Stability:

HTLV-1 Mosaic although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.