

EBV EBNA1 Mosaic Recombinant Protein



RPPB5616

Product Information Protein Information

Product SKU:

RPPB5616

Protein description:

The E.Coli derived recombinant mosaic protein contains the HHV-4 EBNA regions, 1-90, 408-498 amino acids and fused to a 6 aa His Tag at C-terminus and having a molecular weight of 44.2kDa.

Formulation:

10mM PBS pH 7.6 and 10mM NaCl.

Purity:

Protein is >95% pure as determined by SDS-PAGE (coomassie staining).

Stability:

EBV EBNA1 Mosaic protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.