

# EBV EA Recombinant Protein



RPPB5618

---

## Product Information    Protein Information

**Product SKU:**

RPPB5618

**Protein description:**

The E.Coli derived 35.5 kDa recombinant protein contains the HHV-4 Early Antigen Type D, C-terminus regions amino acids 306-390. The protein is fused to a 26kDa GST tag.

**Formulation:**

50mM Tris-HCl, 60mM NaCl, 50% glycerol and 10mM glutathione.

**Purity:**

EBV-Ea protein is >95% pure as determined by SDS-PAGE (coomassie staining).

**Stability:**

EBV-Ea although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.