

# EBV p18 Recombinant Protein



RPPB5619

---

## Product Information    Protein Information

**Product SKU:**

RPPB5619

**Protein description:**

The E.Coli derived recombinant mosaic protein contains the HHV-4 p18 regions, 1-119 amino acids and fused to a GST-Tag at C-terminus.

**Formulation:**

(1mg/ml) containing 25mM Tris-HCl pH 8.0, 1.5M urea and 50% glycerol.

**Purity:**

EBV-p18 protein is >95% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

EBV p18 M although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.