

EBV p18 Recombinant Protein



RPPB5620

Product Information Protein Information

Product SKU:

RPPB5620

Host:

Escherichia Coli.

Protein description:

The E.Coli derived recombinant protein contains the HHV-4 p18 regions, having a molecular mass of 35.4kDa which includes GST tag at N-terminus.

Formulation:

3.8mg/ml in 50mM Tris-HCl, pH 8.0.

Purity:

Protein is >90% pure as determined by 10% PAGE (coomassie staining).

Stability:

EBV p18, GST although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.

Amino Acid Sequence:

ASAGTGALAS SAPSTAVAQS ATPSVSSSIS SLRAATSGAT AAASAAAVD TGSGGGGQPH DTAPRGARKK Q.