

# EBV p18 M Recombinant Protein



RPPB5621

---

## Product Information    Protein Information

**Product SKU:**

RPPB5621

**Host:**

Escherichia Coli.

**Protein description:**

The E.Coli derived recombinant artificial mosaic protein contains two epitopes from C-end associated in one molecule from the HHV-4 p18 regions.

**Formulation:**

50mM Tris-Hcl pH 8.0, 60mM NaCl, 10mM glutathione and 50% glycerol.

**Purity:**

EBV-p18 protein is >95% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

EBV p18 M although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.