

EBV p23 Recombinant Protein



RPPB5622

Product Information Protein Information

Product SKU:

RPPB5622

Host:

Escherichia Coli.

Protein description:

The E.Coli derived recombinant protein contains the HHV-4 p23 regions, 1-162 amino acids and fused to a GST-Tag at C-terminus.

Formulation:

25mM glycine pH-9.6 and 50% glycerol.

Purity:

EBV- p23 protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

EBV-p23 protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.