

# TBEV gE Recombinant Protein



RPPB5625

---

## Product Information   Protein Information

**Product SKU:**

RPPB5625

**Host:**

Escherichia Coli.

**Protein description:**

The E.coli derived recombinant protein contains the Tick-borne Encephalitis Virus glycoprotein E regions, 95-229 amino acids.

**Formulation:**

20mM Tris-MES pH 6.5, 8M urea, 200mM NaCl and 0.05% Tween-20.

**Purity:**

Protein is >90% pure as determined by SDS- PAGE.

**Stability:**

Encephalitis protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.