

TBEV gE middle Recombinant Protein



RPPB5627

Product Information Protein Information

Product SKU:

RPPB5627

Host:

Escherichia Coli.

Protein description:

The E.coli derived recombinant protein contains the Tick-borne Encephalitis Virus glycoprotein E middle regions, 50-250 amino acids.

Appearance:

Sterile Filtered clear solution.

Formulation:

20mM MES pH 6.5, 8M urea, 200mM NaCl and 0.05% Tween-20.

Purity:

Encephalitis protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

Encephalitis protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.