

JEV JEV Recombinant Protein



RPPB5629

Product Information Protein Information

Product SKU:

RPPB5629

Host:

Escherichia Coli.

Protein description:

The E.coli derived recombinant protein contains full length Japanese Encephalitis virus ENV antigen having an Mw of 50kDa. The protein is fused to a 6 histidines tag. GenBank: AHK05344.1

Formulation:

20mM Tris-MES pH 6.5, 8M urea, 200mM NaCl and 0.05% Tween-20.

Purity:

Encephalitis protein is >90% pure as determined by SDS-PAGE (coomassie staining).

Stability:

Encephalitis protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.