

TBEV Core Recombinant Protein



RPPB5630

Product Information Protein Information

Product SKU:

RPPB5630

Protein description:

The E.coli derived recombinant protein contains the Tick-borne Encephalitis Virus core protein epitopes.

Host:

Escherichia Coli.

Formulation:

20mM MES pH 6.5, 8M urea, 200mM NaCl & 0.05% Tween-20.

Purity:

Encephalitis protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

Encephalitis protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.