

# TBEV NS3 Recombinant Protein



RPPB5631

---

## Product Information    Protein Information

**Product SKU:**

RPPB5631

**Protein description:**

The E.coli derived recombinant protein contains the Tick-borne Encephalitis Virus NS3 protein epitopes.

**Host:**

Escherichia Coli.

**Formulation:**

20mM MES pH 6.5, 8M urea, 200mM NaCl and 0.05% Tween-20.

**Purity:**

Encephalitis protein is >95% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

Encephalitis protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.