

# Lassa Capsid Recombinant Protein



RPPB5633

---

## Product Information    Protein Information

**Product SKU:**

RPPB5633

**Host:**

Escherichia Coli.

**Protein description:**

The E.Coli derived Recombinant Lassa Nucleoprotein (LASV strain AV, GenBank accession number AAG41803) containing 420 amino acids of Lassa virus nucleoprotein. The Lassa Capsid protein is fused to a 6xHis tag at C-terminus and purified by proprietary chromatographic technique.

**Appearance:**

Sterile Filtered solution.

**Formulation:**

Lassa Capsid protein solution contains Phosphate buffered saline and 25mM arginine.

**Purity:**

Lassa Capsid is >95% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

Lassa Capsid protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.