

# Rubella E1 Recombinant Protein



RPPB5635

---

## Product Information   Protein Information

**Product SKU:**

RPPB5635

**Protein description:**

The E.Coli derived recombinant protein contains the Rubella Virus E1 regions, 157-176, 374-390, 213-239 amino acids. The protein is fused to a GST-tag.

**Formulation:**

1xpbs.

**Purity:**

Protein is >90% pure as determined by SDS-PAGE.

**Stability:**

Rubella E1 Protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.