

# Rubella E2 Recombinant Protein



RPPB5636

---

## Product Information    Protein Information

**Product SKU:**

RPPB5636

**Protein description:**

The E.Coli derived recombinant protein contains the Rubella Virus E2 regions, 31-105 amino acids.

**Formulation:**

20mM imidazol, 8M urea and 0.3M NaCl.

**Purity:**

Rubella protein is >95% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

Rubella protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.