

Rubella Capsid Recombinant Protein



RPPB5637

Product Information Protein Information

Product SKU:

RPPB5637

Protein description:

The E.Coli derived recombinant protein contains the Rubella Virus Capsid C regions, 1-123 amino acids and fused to a GST-Tag at C-terminus.

Formulation:

25mM Tris-HCl pH 8, 10mM glutathione, 0.2% sarcosyl and 50% glycerol.

Purity:

Rubella protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

Rubella Capsid protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.