

CHIKV E2 Recombinant Protein



RPPB5649

Product Information Protein Information

Product SKU:

RPPB5649

Host:

Escherichia Coli.

Protein description:

Recombinant Chikungunya E2 produced in E.coli having a molecular weight of 38kDa (migrates at 38-40kDa on 10% SDS-PAGE).

Appearance:

Sterile filtered colorless solution.

Formulation:

Sterile Filtered solution containing PBS and 25mM K₂CO₃.

Purity:

Protein is >95% pure as determined by SDS-PAGE.

Stability:

CHIKV E2 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.