

# Zika Envelope F.Length Recombinant Protein



RPPB5651

---

## Product Information    Protein Information

**Product SKU:**

RPPB5651

**Host:**

Escherichia Coli.

**Protein description:**

The E.Coli derived Recombinant Zika Envelope, strain MR766 protein 40-400 a.a. having an Mw of 38kDa. The Zika Envelope full length protein without a C-terminal hydrophobic region. The Zika Envelope full length protein is fused to a 6xHis tag at C-terminus and purified by proprietary chromatographic technique.

**Appearance:**

Sterile Filtered solution.

**Formulation:**

Zika Envelope protein solution contains PBS and 25mM arginine.

**Purity:**

Zika Envelope protein is >95% pure as determined by SDS-PAGE.

**Stability:**

Zika Envelope protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.