

Zika NS1 Recombinant Protein



RPPB5653

Product Information Protein Information

Product SKU:

RPPB5653

Accession:

B5AHE9

Host:

Escherichia Coli.

Protein description:

The E.Coli derived Recombinant Zika NS1 protein (Strain: ZikaSPH2015) having an Mw of 45kDa is derived from the full length Zika NS1 protein. The Zika NS1 protein is fused to a 6xHis tag at C-terminus and purified by proprietary chromatographic technique.

Appearance:

Sterile Filtered solution.

Formulation:

Zika NS1 protein solution contains PBS and 100mM arginine.

Purity:

Zika NS1 protein is >95% pure as determined by SDS-PAGE.

Stability:

Zika NS1 protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.