

Zika NS1 Recombinant Protein



RPPB5654

Product Information Protein Information

Product SKU:

RPPB5654

Host:

Human Embryonic
Kidney 293 cells.

Protein description:

The HEK derived recombinant Zika NS1 protein strain MR-766 (Prototype, Uganda, 1947, GenBank accession number AY632535). The Zika NS1 protein is fused to a 6xHis tag at C-terminus and purified by proprietary chromatographic technique.

Appearance:

Sterile Filtered solution.

Formulation:

Zika NS1 protein contains Dulbecco's Phosphate buffered saline, pH 7.4.

Purity:

Zika NS1 protein is >90% pure as determined by SDS-PAGE.

Stability:

Zika NS1 protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.