

Zika NS1 Recombinant Protein



RPPB5655

Product Information Protein Information

Product SKU:

RPPB5655

Host:

Sf9, Baculovirus cells.

Protein description:

Sf9 derived recombinant Zika NS1 protein. The Zika NS1 protein is fused to a 6xHis tag and purified by proprietary chromatographic technique.

Appearance:

Sterile Filtered solution.

Formulation:

Zika NS1 protein contains Phosphate buffered saline, pH 7.4 and 0.09% NaN₃.

Purity:

Zika NS1 protein is >85% pure as determined by SDS-PAGE.

Stability:

Zika NS1 protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.