

# Zika Envelope Recombinant Protein



RPPB5656

---

## Product Information    Protein Information

**Product SKU:**

RPPB5656

**Host:**

Sf9, Baculovirus cells.

**Protein description:**

Sf9 derived recombinant Zika Envelope protein, having an Mw of approximately 50kDa. The Zika Envelope protein is fused to a 6xHis tag and purified by proprietary chromatographic technique.

**Appearance:**

Sterile Filtered solution.

**Formulation:**

Zika Envelope protein contains Phosphate buffered saline, pH 7.4 and 0.09% NaN<sub>3</sub>.

**Purity:**

Zika Envelope protein is >80% pure as determined by SDS-PAGE.

**Stability:**

Zika Envelope protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.