

# Parvovirus B19 VLP VP1/VP2 Recombinant Protein



RPPB5660

---

## Product Information    Protein Information

**Product SKU:**

RPPB5660

**Host:**

Sf9 insect cells.

**Protein description:**

Recombinant Parvovirus B19 VLP VP1/VP2 Co-Capsid produced in SF9 containing B19 capsid proteins VP1 (81kDa) and VP2 (61kDa). Parvovirus B19 VLP VP1/VP2 is purified by proprietary chromatographic techniques.

**Appearance:**

Sterile Filtered clear solution.

**Formulation:**

Parvovirus B19 VLP VP1/VP2 protein solution is supplied in 16mM Sodium phosphate pH-7.4, 120mM NaCl and 20% glycerol.

**Stability:**

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.