

Parvovirus B19 VLP VP1 Recombinant Protein



RPPB5661

Product Information Protein Information

Product SKU:

RPPB5661

Host:

e.coli.

Protein description:

Recombinant Parvovirus B19 VLP VP1 produced in e.coli contains 227 amino acids from the N-terminus of VP1 protein sequence. Parvovirus B19 VLP VP1/VP2 is purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered clear solution.

Formulation:

Parvovirus B19 VLP VP1 protein solution is supplied in PBS pH 7.4 and 0.02% Sodium azide.

Purity:

Greater than 90.0% as determined by 12% PAGE (coomassie staining).

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.