

# H5N1 Indonesia Recombinant Protein



RPPB5665

---

## Product Information    Protein Information

**Product SKU:**

RPPB5665

**Host:**

Baculovirus Insect  
Cells.

**Protein description:**

Recombinant Full-Length H5N1 A/Indonesia/05/2005 is glycosylated with N-linked sugars, produced using baculovirus vectors in insect cells.

**Appearance:**

Sterile Filtered colorless solution.

**Formulation:**

The Recombinant H5N1 A/Indonesia/05/2005 solution contains 10mM Sodium phosphate, pH 7.2, 150mM NaCl and 0.005% Tween-20.

**Purity:**

Greater than 90.0% as determined by SDS-PAGE.

**Stability:**

H5N1 A/Indonesia/05/2005 Recombinant should be stored at 4°C. Do NOT Freeze!